

# Rubber Fab Technologies Group's **GAUGE GUARD ISOLATOR HYGIENIC SEALS**

Protect Critical Instrumentation  
With No Performance Loss



The Rubber Fab Gauge Guard Isolator seals **protect expensive diaphragm and instruments from damage without affecting instrument performance**. By combining a quality Rubber Fab hygienic seal with a membrane, the Gauge Guard Isolator hygienic seal **isolates hot or caustic saline solution from instrumentation**. A necessary and **cost effective** step in **assuring long life and accurate results from your instruments and gauges**. As some manufacturers will suggest the use of a gauge guard seal, others note that the **failure to use a gauge guard seal can void most instrument and gauge warranties**.



## All the Rubber Fab Benefits in a Gauge Guard Hygienic Seal

- Seal membrane **protects instrumentation from caustic saline solutions**
- **Will not interfere with gauge operation or accuracy**
- Works with most industry standard instruments
- Helps **extend the life** of gauges
- Available in:  
**FKM Fluoroelastomer**  
**Platinum Silicone**  
**EPDM**
- 1/2" and 3/4" – .015 membrane thickness
- 1", 1-1/2" and 2" – .010 membrane thickness
- Metal Detectable elastomers available

## All Rubber Fab Hygienic Seals Meet the Stringent Standards For Purity

- U.S. Pharmacopeia Class VI Certification
- Cytotoxicity Criteria
- Title 21 CFR 177.2600
- USDA and 3-A Sanitary Standards
- Current Good Manufacturing Practices (CGMP)
- ASME-BPE Standards
- Animal Derived Ingredient Free

© ADI Free is a registered trademark of Rubber Fab Technologies Group.

Rubber Fab Technologies Group's Metal Detectable products are patent pending.




Rubber Fab Technologies Group  
Limecrest Road

P.O. Box 626 • Andover, NJ 07821

phone: 973-579-2959 fax: 973-579-7275

[www.rubberfab.com](http://www.rubberfab.com)




### FKM Fluoroelastomer

1/2"	42GGI-SFY-050
3/4"	42GGI-SFY-075
1"	40GGI-SFY-100
1-1/2"	40GGI-SFY-150
2"	40GGI-SFY-200



### Platinum Silicone

1/2"	42GGI-PX-050
3/4"	42GGI-PX-075
1"	40GGI-PX-100
1-1/2"	40GGI-PX-150
2"	40GGI-PX-200



### EPDM

1/2"	42GGI-E-050
3/4"	42GGI-E-075
1"	40GGI-E-100
1-1/2"	40GGI-E-150
2"	40GGI-E-200

## New From Rubber Fab Technologies Group

**DESIGN**  
validation  
production  
**QUALITY ASSURANCE**

Improving Product Quality and Consistency

**RUBBER FAB**

**Design...Validation...Production  
...Quality Assurance**  
*Improving Product Quality and Consistency.*

Rubber Fab Technologies Group is a member of:

Distributed by:



International Association of  
Food Industry Suppliers



Distributors Serving  
Sanitary Processing

